

Scan Body Reference Guide

Digitizing Cases with Scan Bodies



**Your Connection
for Every Connection**

azure[™]
Multi-Platform Prosthetic Solutions

Scan Bodies for Azure Multi-Platform Prosthetic Solutions

- Digitize restorative cases utilizing Azure MPS components in RealGUIDE™, 3Shape, Exocad™, and Dental Wings
- Utilize scan bodies from multiple manufacturers for 20 implant systems
- To access, update your libraries from <https://azuredental.com/resources/digital-library.html>

Scan Body Library Availability

Implant System	Azure	Atlantis™	Elos	3Shape*	Straumann® Mono	Straumann® Metal
ASTRA TECH IMPLANT SYSTEM™ EV	•	•	•	•		
ASTRA TECH OSSEOSPEED™	•	•	•	•		
AZURE TAPERED BONE LEVEL	•	•	•	•	•	•
BIOHORIZONS® INTERNAL	•	•	•	•		
DENTIUM SUPERLINE™ & IMPLANTIUM®	•					
DENTSPLY ANKYLOS® C/X	•	•	•	•		
MEGAGEN ANYRIDGE®	•					
MIS® SEVEN	•	•				
NEODENT® GRAND MORSE®	•		•			
NOBELACTIVE® / REPLACE®	•	•	•	•		
NOBEL BRÅNEMARK®	•		•	•		
NOBELREPLACE SELECT™	•	•	•	•		
OSSTEM® TS / HIOSSEN® ET	•	•	•	•		
STRAUMANN® BLX	•					•
STRAUMANN® BONE LEVEL	•	•	•	•	•	•
STRAUMANN® TISSUE LEVEL	•		•	•	•	•
ZIMVIE CERTAIN®	•	•	•	•		
ZIMVIE EXTERNAL HEX	•		•	•		
ZIMVIE TSV®	•	•	•	•		
MULTI-UNIT®	•		•			

*only for use in 3Shape systems

Azure offers scan bodies for all 30+ implant systems. However, scan bodies from other third-party systems can also be utilized.

Neodent®	NT-Trading	Medentika®	TruAbutment™	Preat™ PEEK Blue	Preat™ Titanium	CEREC® Omnicam	CEREC® Bluecam	CEREC® ScanPost

			
			
			.		.			

.		.						

	3.0 mm(P)	3.6 mm(P)	4.2 mm(P)	4.8 mm(P)	3.6 mm(P)
DAEV-CONICAL compatible with Astra Tech Implant System™ EV	EV	EV	EV	EV	UniAbutment
AZURE	ARDAEV-ISB3	ARDAEV-ISB36	ARDAEV-ISB42	ARDAEV-ISB48	ARDAEVMU-ISBNH36
ATLANTIS	IO-P-01 (35243)	IO-P-02 (35244)	IO-P-03 (35245)	IO-P-04 (35246)	-
ELOS	IO 3B-A	IO 3B-B	IO 3B-C	IO 3B-D	IO 3D-A
3SHAPE*	-	3S02-DEP	3S02-DEY	3S02-DEB	-
NT-TRADING	SEV 9.S3D3.000	SEV 9.S3D3.600	SEV 9.S3D4.200	SEV 9.S3D4.800	-
MEDENTIKA	EV1400	EV1410	EV1420	EV1430	-
TRUABUTMENT	AE30-SB	AE36-SB	AE42-SB	AE48-SB	-
PREAT TITANIUM	9006445	9006545	9006645	9006745	-
CEREC OMNICAM	-	6431311	6431329	6431329	-
CEREC BLUECAM	-	6431295	6431303	6431303	-
CEREC SCANPOST	-	6586361	6586379	6586387	-

*only for use in 3Shape systems

	3.0 mm(P)	3.5/4.0 mm(P)	4.5/5.0 mm(P)	3.6 mm(P)	3.6 mm(P)
DATX-INTERNAL compatible with Astra Tech Osseospeed™	NP	RP	WP	UniAbutment 20°	UniAbutment 40°
AZURE	ARDATX-ISB3	ARDATX-ISB35	ARDATX-ISB45	ARDATX20-ISBNH36	ARDATX45-ISBNH36
ATLANTIS	IO-C-01 (35215)	IO-C-02 (35216)	IO-C-03 (35217)	-	-
ELOS	IO 3A-A	IO 3A-B	IO 3A-C	IO 3D-B	IO 3D-C
3SHAPE*	3S02-DAX	3S02-DAS	3S02-DAL	-	-
NT-TRADING	-	S9.S3D3.540	S9.S3D4.550	-	-
MEDENTIKA	S1430	S1400	S1420	-	-
TRUABUTMENT	AS3-SB	AS4-SB	AS5-SB	-	-
CEREC OMNICAM	-	6431329	6431329	-	-
CEREC BLUECAM	-	6431303	6431303	-	-
CEREC SCANPOST	-	6431055	6431063	-	-

*only for use in 3Shape systems

	3.0 mm(P)	3.5 mm(P)	4.5 mm(P)	5.7 mm(P)
BHI-INTERNAL compatible with BioHorizons® Internal	-	NP	RP	WP
AZURE	ARBHI-ISB3	ARBHI-ISB35	ARBHI-ISB45	ARBHI-ISB57
ATLANTIS	-	IO-B-01 (35212)	IO-B-02 (35213)	IO-B-03 (35214)
ELOS	-	IO 5A-A	IO 5A-B	IO 5A-C
3SHAPE*	-	3S02-BIY	3S02-BIG	3S02-BIB
NT-TRADING	BIO9.S3D3.000	BIO9.S3D3.500	BIO9.S3D4.500	BIO9.S3D5.700
MEDENTIKA	-	R1400	R1410	R1420
TRUABUTMENT	BHN-SB	BH35-SB	BH45-SB	BH57-SB
CEREC OMNICAM	-	6431329	6431329	6431329
CEREC BLUECAM	-	6431303	6431303	6431303
CEREC SCANPOST	-	6532886	6532944	6536234

*only for use in 3Shape systems

	3.6 mm(P)	3.8 mm(P)	4.0 mm(P)	4.5 mm(P)	5.0 mm(P)
DT-CONICAL compatible with Dentium Superline™ & Implantium®	-	RP	-	-	-
AZURE			ARDT-ISBU		
MEDENTIKA			DT1400L		
TRUABUTMENT			DSU-SB		

DA-INTERNAL compatible with Dentsply Ankylos® C/X	
AZURE	ARDA-ISB35
ATLANTIS	IO-L-01
ELOS	IO 8B-A
3SHAPE*	3S02-DYR
NT-TRADING	Y9.S3D3.50X
MEDENTIKA	Y1400X-L
CEREC OMNICAM	6431311
CEREC BLUECAM	6431295
CEREC SCANPOST	6586569

*only for use in 3Shape systems

MA-CONICAL compatible with Megagen AnyRidge®	
AZURE	ARMA-ISBU
NT-TRADING	MEG9.S3D3.555
MEDENTIKA	MG1400-L
TRUABUTMENT	MRD-SB
PREAT PEEK BLUE	8004440
PREAT TITANIUM	9004445

3.3 mm(P) 3.75/4.2 mm(P) 5.0/6.0 mm(P)

MSH-INTERNAL compatible with Mis® Seven	3.3 mm(P)	3.75/4.2 mm(P)	5.0/6.0 mm(P)
	NP	SP	WP
AZURE	ARMSH-ISB33	ARMSH-ISB34	ARMSH-ISB56
ATLANTIS	IO-B-04 (35273)	IO-B-01 (35212)	IO-B-02 (35213)
NT-TRADING	MIS9.S3D3.300	MIS9.S3D3.742	MIS9.S3D5.000
MEDENTIKA	R1400	R1410	R1420

NG-CONICAL compatible with Neodent® Grand Morse®	3.3 mm(P)
	GM
AZURE	ARNG-ISB3
ELOS	IO 13A-A
TRUABUTMENT	NDG-SB
PREAT TITANIUM	9007145

3.0 mm(P) 3.5 mm(P) 4.3/5.0 mm(P) 5.5 mm(P)

NAC-CONICAL compatible with NobelActive® / Replace®	3.0 mm(P)	3.5 mm(P)	4.3/5.0 mm(P)	5.5 mm(P)
	-	NP	RP	WP
AZURE	ARNAC-ISB3	ARNAC-ISB35	ARNAC-ISB43	ARNAC-ISB55
ATLANTIS	IO-K-03 (35266)	IO-K-01 (35233)	IO-K-02 (35234)	IO-K-04 (35727)
ELOS	IO 2B-D SA	IO 2B-A SA	IO 2B-B SA	IO 2B-C SA
3SHAPE*	3S02-NA3	3S02-NAN	3S02-NAR	3S02-NAW
NT-TRADING	-	F9.S3D3.500	F9.S3D4.350	-
MEDENTIKA	F1420	F1400	F1410	-
TRUABUTMENT	NA3-SB	NAN-SB	NAR-SB	NAW-SB
PREAT PEEK BLUE	-	8000840	8000940	-
PREAT TITANIUM	9004745	9000845	9000945	-
CEREC OMNICAM	-	6431329	6431329	-
CEREC BLUECAM	-	6431303	6431303	-
CEREC SCANPOST	-	6431279	6431287	-

*only for use in 3Shape systems

	3.5 mm(P)	4.1 mm(P)	5.1 mm(P)
NB-EXTERNAL compatible with Nobel Brånemark®	NP	RP	WP
AZURE	ARNB-ISB35	ARNB-ISB41	ARNB-ISB51
ELOS	IO 6A-A	IO 6A-B	IO 6A-C
3SHAPE*	3S02-NBN	3S02-NBR	3S02-NBW
NT-TRADING	K9.S3D3.500	K9.S3D4.100	K9.S3D5.100
MEDENTIKA	K1400	K1410	K1420
CEREC OMNICAM	6431329	6431329	-
CEREC BLUECAM	6431303	6431303	-
CEREC SCANPOST	6431006	6431022	-

*only for use in 3Shape systems

	3.5 mm(P)	4.3 mm(P)	5.0 mm(P)	6.0 mm(P)
NRT-TRILOBE compatible with NobelReplace Select™	NP	RP	WP	6.0
AZURE	ARNRT-ISB35	ARNRT-ISB43	ARNRT-ISB5	ARNRT-ISB6
ATLANTIS	IO-D-01 (35219)	IO-D-02 (35220)	IO-D-03 (35221)	IO-D-04 (35222)
ELOS	IO 2A-A	IO 2A-B	IO 2A-C	IO 2A-D
3SHAPE*	3S02-NRN	3S02-NRR	3S02-NAR	3S02-NRW
NT-TRADING	E9.S3D3.500	E9.S3D4.300	E9.S3D5.000	E9.S3D6.000
MEDENTIKA	E1400	E1410	E1420	E1430
TRUABUTMENT	NRN-SB	NRR-SB	NRW-SB	NR6-SB
PREAT PEEK BLUE	8000440	8000540	8000640	-
PREAT TITANIUM	9000445	9000545	-	-
CEREC OMNICAM	6431329	6431329	6431329	6431329
CEREC BLUECAM	6431303	6431303	6431303	6431303
CEREC SCANPOST	6430933	6430941	6430958	6430982

*only for use in 3Shape systems

	Mini	Regular
OH-CONICAL compatible with Osstem® TS / HIOSSEN® ET	-	STD
AZURE	AROH-ISB3	AROH-ISB4
ATLANTIS	IO-R-01 (35270)	IO-R-02 (35271)
ELOS	IO 12A-A	IO 12A-B
3SHAPE*	3S02-OTM	3S02-OTR
NT-TRADING	HIO9.S3D3.500	HIO9.S3D4.060
MEDENTIKA	OT1400-L	OT1410-L
TRUABUTMENT	OTM-SB	OTR-SB
PREAT PEEK BLUE	8005040	8005140
PREAT TITANIUM	9005045	9005145
CEREC OMNICAM	6431329	6431329
CEREC BLUECAM	6431303	6431303
CEREC SCANPOST	6534197	6536846

*only for use in 3Shape systems

SBLX-CONICAL compatible with Straumann® BLX	RB/WB
AZURE	ARSBLX-ISB35
STRAUMANN® METAL	065.0103
MEDENTIKA	LX1400
TRUABUTMENT	SXR-SB
PREAT TITANIUM	9007045

	3.3 mm(P)	4.1/4.8 mm(P)
SBL-CONICAL compatible with Straumann® Bone Level & Azure Tapered Bone Level	NC	RC
AZURE	ARSBL-ISB33	ARSBL-ISB41
ATLANTIS	IO-H-01 (35228)	IO-H-02 (35229)
ELOS	IO 4A-A	IO 4A-B
3SHAPE*	3S02-SBN	3S02-SBR
STRAUMANN® MONO	025.2915	025.4915
STRAUMANN® METAL	025.0078	025.0079
NT-TRADING	L9.S3D3.30	L9.S3D4.148
MEDENTIKA	L1400	L1410
TRUABUTMENT	SBN-SB	SBR-SB
PREAT PEEK BLUE	8001540	8001440
PREAT TITANIUM	9001545	9001445
CEREC OMNICAM	6431329	6431329
CEREC BLUECAM	6431303	6431303
CEREC SCANPOST	6431246	6431253

*only for use in 3Shape systems

	3.5 mm(P)	4.1/4.8 mm(P)	4.8/6.5 mm(P)
STL-TISSUE compatible with Straumann® Tissue Level	NNC	NC	WN
AZURE	ARSTL-ISB35	ARSTL-ISB41	ARSTL-ISB48
ELOS	IO 4C-A	IO 4B-A	IO 4B-B
3SHAPE*	3S02-SSN	3S02-SSR	3S02-SSW
STRAUMANN® MONO	048.173	048.168	048.169
STRAUMANN® METAL	048.384	048.385	-
NT-TRADING	N9.S3D3.500	N9.S3D4.800	N9.S3D6.500
MEDENTIKA	N1430	N1410	N1420
TRUABUTMENT	-	STR-SB	STW-SB
PREAT PEEK BLUE	-	8001240	-
PREAT TITANIUM	-	9001245	-
CEREC OMNICAM	-	6431329	6431329
CEREC BLUECAM	-	6431303	6431303
CEREC SCANPOST	-	6431170	6431196

*only for use in 3Shape systems

	3.4 mm(P)	4.1 mm(P)	5.0 mm(P)
ZC-CERTAIN® compatible with ZimVie CERTAIN®	NP	RP	WP
AZURE	ARZC-ISB34	ARZC-ISB41	ARZC-ISB5
ATLANTIS	IO-F-01 (35225)	IO-F-02 (35226)	IO-F-02 (35226)
ELOS	IO 7B-A	IO 7B-B	IO 7B-B
3SHAPE*	3S02-BCP	3S02-BCB	3S02-BCY
NT-TRADING	H9.S3D3.400	H9.S3D4.150	H9.S3D4.150
MEDENTIKA	H1400	H1410	H1410
TRUABUTMENT	BC34-SB	BC41-SB	BC50-SB
PREAT PEEK BLUE	8001640	8001740	-
PREAT TITANIUM	9001645	9001745	-
CEREC OMNICAM	6431311	6431329	6431329
CEREC BLUECAM	6431295	6431303	6431303
CEREC SCANPOST	6431212	6431220	6431238

*only for use in 3Shape systems

	3.4 mm(P)	4.1 mm(P)	5.0 mm(P)
ZEX-EXTERNAL compatible with ZimVie EXTERNAL HEX	NP	RP	WP
AZURE	ARZEX-ISB34	ARZEX-ISB41	ARZEX-ISB5
ELOS	IO 7A-A	IO 7A-B	IO 7A-B
3SHAPE*	3S02-BEP	3S02-BEB/3S02-BE4	3S02-BEY
NT-TRADING	I 9.S3D4.150	I 9.S3D4.150	I 9.S3D5.000
MEDENTIKA	I1400	I1410	I1410
CEREC OMNICAM	6431329	6431329	6431329
CEREC BLUECAM	6431303	6431303	6431303
CEREC SCANPOST	6431089	6431105	6431113

*only for use in 3Shape systems

	3.5 mm(P)	4.5 mm(P)	5.7 mm(P)
ZTSV-INTERNAL compatible with ZimVie TSV®	NP	RP	WP
AZURE	ARZTSV-ISB35	ARZTSV-ISB45	ARZTSV-ISB57
ATLANTIS	IO-B-01 (35212)	IO-B-02 (35213)	IO-B-03 (35214)
ELOS	IO 5A-A	IO 5A-B	IO 5A-C
3SHAPE*	3S02-ZTG	3S02-ZTP	3S02-ZTY
NT-TRADING	R9.S3D3.500	R9.S3D4.500	R9.S3D5.700
MEDENTIKA	R1400	R1410	R1420
TRUABUTMENT	ZV3-SB	ZV4-SB	ZV5-SB
PREAT PEEK BLUE	8000140	8000240	8000340
PREAT TITANIUM	9000145	9000245	9000345
CEREC OMNICAM	6431329	6431329	6431329
CEREC BLUECAM	6431303	6431303	6431303
CEREC SCANPOST	6431139	6431147	6431238

*only for use in 3Shape systems

	4.8 mm(P)	6.5 mm(P)
Multi-Unit®	RP	WP
AZURE	ARMU-ISB48	ARMU-ISB65
ELOS	IO 2C-A	IO 2C-B
NEODENT	108.218	-
MEDENTIKA	8500	-

azure™ | A ZimVie Brand

For more information, visit azuredental.com

Dental Wings and DWOS® are registered trademarks of Dental Wings Inc, Exocad™ is a trademark of exocad GmbH, ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XIVE®, FRIADENT®, ASTRATECH implant system™ EV and UniAbutment®, Cerec® and Atlantis are trademarks and/or registered trademarks of DENTSPLY SIRONA Implants BioHorizons® is a registered trademark of BioHorizons, Inc. Bränemark System®, NobelReplace® and Replace Select™, NobelActive™, NobelReplace® NobelSpeedy™ and Multi-Unit® are trademarks and/or registered trademarks of Nobel Biocare AB Neodent®, Grand Morse®, Medentika® and Straumann® are registered trademarks of Straumann Holding AG, Switzerland MIS® is a registered trademark of MIS Implant Technologies, Ltd. Tapered Screw-Vent® and SwissPlus® are registered trademarks of Zimmer Dental, Inc. USA Certain® is a registered trademark of Biomet 3i Camlog® and Conelog® are registered trademarks of CAMLOG Biotechnologies AG Megagen® is a registered trademark of MEGAGEN IMPLANT CO. LTD Osstem® and Hiossen® are registered trademarks of OSSTEM IMPLANT CO. LTD Dentium® is a registered trademark of Dentium Co., Ltd. Preat is a registered Trademark of Preat Corporation Elos is a registered trademark of Elos Medtech AB 3Shape is a registered trademark of 3Shape A/S NT Trading is a registered trademark of NT-Trading GmbH & Co. KG TruAbutment is a registered trademark of TruAbutment Inc

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. Azure MPS abutments, components and tools are manufactured by Terrats Medical S.L. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV2636 REV A 03/25 ©2025 ZimVie Inc. All rights reserved.

